

REJUBEAU HA STYLISH PLUS

REJUBEAU: BEAUTY IS ON

REJUBEAU HA STYLISH PLUS with cutting edge Technology: MSCT



About REJUBEAU HA STYLISH PLUS

REJUBEAU HA STYLISH PLUS is not only safe, stylish, & long lasting, but also has excellent molding ability based on MSCT (Multi Staged Cross-linking Technology)



Safe



Long Lasting




Excellent Molding Ability








Stylish

About the Technology of REJUBEAU HA STYLISH PLUS




MSCT Multi Staged Cross-linking Technology




A. Comparison of Cross-linking Agent (BDDE vs DVS)

 Cross-linking Technology	CONTENT	BDDE	DVS
	Chemical Name	1,4-Butanediol diglycidyl ether	Divinyl Sulfone
	Structure		
	Molecular Weight (Dalton)	202.25	118.15
	Molecule Length (Ratio)	3 (in case)	1 (then approx. 1/3)
	The number of atoms forming a bridge between HA molecules	14 units (Possibility of surplus atoms ↑)	5 units
	Manufacturing Characteristics	Easy Process Loose Cross-linking	Difficult Process Dense Cross-linking
	Cross-linking Degree	Due to incomplete Cross-linking pendent phenomenon occurs	Complete Cross-linking
	Cross-link		

B. Property comparison of HA Filler (Biphasic vs Monophasic)

 Property Composition	CONTENT	Diagram	Properties	Feature
	Bi-Phasic		HA : Exist in 2 physical properties (Gel +Liquid)	<ul style="list-style-type: none"> - High gel fluidity - Low viscosity ; easy shooting - May flow down and disperse
Mono-Phasic		HA : 1 Cross-linked physical property (Gel)	<ul style="list-style-type: none"> - Excellent gel cohesion and viscoelasticity - Requires skill for shooting - Neither flow down nor disperse 	

Longevity comparison graph




Variations to conventional products that claim to have produced a long-lasting HA filler:

Most companies claim to have successfully manufactured a long-lasting HA filler by increasing the Cross-linking rate. However, due to limitations in the Cross-linking technology, these products instead have increased viscosity and particle size, which does not actually contribute to a life span of HA filler. Moreover, products that have only raised viscosity fail to demonstrate good moldability.

REJUBEAU HA STYLISH PLUS Test of Viscosity & Molding

The results of measuring the changes of our products and our competitor's products, after a time elapse of four minutes. As shown below, our competitor's products demonstrated observable change while our products proved to maintain its shape.



without altering its chemical properties




Elastic and glossy
REJUBEAU HA STYLISH PLUS Outstanding shape retention

Simple particle enlargement due to limit in technological skills
Reduces longevity of the product


Safe & Removable

Results of measuring the sliding down of our products and other company's products. After 4 minutes elapse as shown below, while third-party products slammed down a lot, our products can be seen to keep shapes on an ongoing basis.




Superior hydrophilic capacity




RUJUBEAU HA STYLISH PLUS provides an excellent volume effect due to its high hydrophilic capacity and stable molecular structure.



Ergonomic Design



Characteristics of REJUBEAU HA STYLISH PLUS

Product	HA Concentration	Volume	Needle	Quality	Viscosity	Molding	Injecting Area
 REJUBEAU HA STYLISH PLUS 300	20mg/mL	1mL X 1 Syringe 0.3% Lidocaine	25G 13mm UTW Needle Cannula 23G50	★★★★	★★★★	★★★★	Upper Dermis
 REJUBEAU HA STYLISH PLUS 200	20mg/mL	1mL X 1 Syringe 0.3% Lidocaine	27G 13mm UTW Needle Cannula 23G50	★★★★	★★★☆☆	★★★☆☆	Lower Dermis, Upper sub-Q
 REJUBEAU HA STYLISH PLUS 100	20mg/mL	1mL X 1 Syringe 0.3% Lidocaine	30G 13mm UTW Needle	★★★★	★★★☆☆	★★★☆☆	Lower sub-Q

*UTW: Ultra Thin Wall



Application of REJUBEAU HA STYLISH PLUS



Forehead

REJUBEAU HA STYLISH PLUS 300

Deep Wrinkles

REJUBEAU HA STYLISH PLUS 200

Fine Wrinkles

REJUBEAU HA STYLISH PLUS 100

Nose

REJUBEAU HA STYLISH PLUS 300

Cheekbone

REJUBEAU HA STYLISH PLUS 200, 300

Nasolabial

REJUBEAU HA STYLISH PLUS 300

Lip

REJUBEAU HA STYLISH PLUS 100, 200

Short Chin

REJUBEAU HA STYLISH PLUS 300



TRACLE INC



801, Samyoong Tower, 1090, Gyeongui-ro, Paju-si,
Gyeonggi-do 10908, Republic of Korea
e. tmicmo@gmail.com / t. +82 70 4696 5346 /
f. +82 31 943 8183 / www.tracle.kr / Instagram: [tracle_korea](https://www.instagram.com/tracle_korea)